




Radial Hinged Ulnar Deviation Splint

To reduce ulnar drift at the MP joints caused by rheumatoid arthritis, neurological impairment or trauma. Use following MP arthroplasty to maintain alignment after removal of post-operative splint. The Radial Hinged Ulnar Deviation Splint can be used to align strongly deviated fingers.




Open finger straps to allow for easy donning of the splint.



Slide hand and fingers into the Radial Hinged splint. Contour the plastic covered stay against the palm to secure the splint on the hand.



Adjust the finger strap as needed to bring the fingers into alignment. Once adjusted, it may be possible to slide the fingers in and out of the splint without reopening the finger strap.



Splint shown on right hand.



Stevensville, MD 21666 U.S.A.

Made in the U.S.A.

ADJUSTMENTS



The RHUD splint is made from 1/16" Kydex. Adjustments should be made by a health care professional. **CAUTION:** Kydex is a high temperature plastic and cannot be placed directly on the skin when hot and pliable. Place stockinette or other covering on hand when contouring warm material. To custom fit the RHUD, use a heat gun to warm the ulnar side until pliable. Cover the wearers hand and contour the splint around the ulnar side of the hand. Trim warm material with scissors if needed. Padding may be added as needed.

TO ORDER

Measure the width across the 2nd to the 5th MP joints.

	Width at MP's	Right	Left
Small	2 3/4" - 3"	P2003-R2	P2003-L2
Medium	3" - 3 3/8"	P2003-R3	P2003-L3
Large	3 3/8" - 3 5/8"	P2003-R4	P2003-L4

Replacement strap available in packages of 2

Order Toll Free **888.378.7763** (U.S only)

Local/Int. 410.604.6393

24 hour Fax 410.604.6398

Monday - Friday 8:30 - 5:30 pm EST

CARE INSTRUCTIONS



Wash with mild soap and water. Strap can be hand washed and air dried.

Material: Kydex, Nylon Strap, cotton and wool padding.

©2005, 3-Point Products, Inc.

Authorized European Representative
Emergo Europe
Molenstraat 15
2513 BH The Hague
The Netherlands



www.3pointproducts.com